Safety and Effectiveness is compromised with more than one use of each of the contents of this kit.

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Report any serious incident in relation to this refill kit to the manufacturer. In the European Union, report any serious incident to the manufacturer and the competent authority of the Member State.

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AN 1004
AN 7514
AN 2018

Ethylene Oxide Cartridges

For use in the following Andersen Sterilizers Ethylene Oxide Gas Sterilization Systems:

Anprolene AN74i/ix/j
Anprolene AN75i/j
EOGas 4
EOGas 3

Manufactured by:
Andersen Sterilizers, Inc.
Health Science Park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.

Distributed in the US by:
Andersen Sterilizers, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 U.S.A.

H.W. Andersen Products of California, Inc.
3151 Caroline Drive
Haw River, NC 27258 U.S.A.

Authorized UK Representative:
H.W. Andersen Products, Ltd.
808 Fowler Road
Clacton-on-Sea • Essex CO15
4AA U.K.

Authorized EU Representative:
H. W. ANDERSEN S.A.R.L.
12 Place Saint-Hubert • 59000 •
Lille, France

For all other inquiries worldwide:

Andersen Sterilizers, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 USA

ECN	Revision	Effective Date
2023032	0	2024-08-16

Cartridges

AN1004

AN7514

AN2018

AN 1004 EOGAS 4®

CONTENTS: 14 Sterilization Bags, 14 AN1036 Dosimeter®, 14 Humidichip®, 14 Cartridges with ampoules

AN7514 ANPROLENE®

CONTENTS: 14 Sterilization Bags, 14 AN87 Dosimeter®, 14 Cartridges with ampoules

AN2018 EOGAS® 3

CONTENTS: 14 Cartridges with ampoules

Content of Cartridges: each ampoule contains 0.64 av. oz. (17.6 g)

Ethylene Oxide

Active ingredient: Ethylene oxide 97.0%

Inert ingredient: 3.0%

Total: 100%

For use in the following Andersen Sterilizers Ethylene Oxide Gas Sterilization Systems:

Anprolene AN74i/ix/j

Anprolene AN75i/j

EOGas 4

EOGas 3

Manufactured by:

ANDERSEN STERILIZERS, INC.

3154 Caroline Drive • Haw River, NC 27258 USA

Distributed by:

ANDERSEN STERILIZERS, INC.

3202 Caroline Drive • Haw River, NC 27258 USA

1-800-523-1276 • 336-376-3000

Authorized UK Representative:

H. W. ANDERSEN PRODUCTS, LTD.

808 Fowler Clacton-On-Sea • Essex CO15 4AA UK

44-1255-428328

Authorized EU Representative:
H. W. ANDERSEN S.A.R.L.
12 Place Saint-Hubert • 59000 • Lille, France

EPA Registration No. 69340-7
EPA Establishment No. 69340-NC-001

Keep Out of Reach of Children

DANGER

PRECAUTIONARY STATEMENTS

Causes eye and skin burns. May be fatal if inhaled. Do not breathe vapor. Do not swallow. Do not get in eyes, on skin, or on clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Wear appropriate protective clothing.

USER SAFETY RECOMMENDATIONS

Users should remove clothing/PPE immediately if pesticide gets inside. Then wash thoroughly and put on clean clothing. Users should remove PPE immediately after handling this product. Wash the outside of gloves before removing. As soon as possible, wash thoroughly and change into clean clothing.

OTHER POSSIBLE DELAYED HEALTH EFFECTS

Cancer and reproductive hazard. May cause nervous system damage, cataracts, adverse reproductive effects, chromosomal and mutagenic changes. May cause irritation of respiratory tract, chest tightness, headache, nausea, vomiting, diarrhea, light headed feeling, dizziness, weakness, drowsiness, cyanosis, loss of coordination, coma, delayed lung injury (fluid in lungs), immediate or delayed skin irritation and blisters, allergic skin.

PEL: 1 PPM-TWA Ethylene Oxide (OSHA-29CFR 1910.1047)
EL: 5 PPM-excursion limit, 15 minutes
Store and use with adequate ventilation in accordance with 29 CFR1910.1047

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FIRST AID

In all cases of exposure, get medical attention immediately. Take person to a doctor or emergency treatment facility at once. Have the product container or label with you when calling the Emergency Contact number or doctor or going for treatment.

IF INHALED: Move exposed person to fresh air. Keep warm. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth method. If breathing is at all labored, give oxygen. Call a doctor even if no symptoms are present for further treatment advice. Keep under medical observation. Symptoms may be delayed.

IF IN EYES: Hold eyelids open and flush eyes with a steady, gentle stream of water for at least 15 to 20 minutes. Remove contact lenses, if present, after the first 5 minutes and then continue rinsing the eyes. Get immediate medical treatment.

IF ON SKIN: Immediately wash skin for 15-20 minutes with plenty of water while removing contaminated clothing and shoes. Call Emergency Contact number or a doctor for treatment advice. Aerate, wash, or clean contaminated clothing and discard leather goods.

IF SWALLOWED: Call the Emergency Contact number or a doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by an Emergency Contact responder or doctor. Do not give anything to an unconscious person.

NOTE TO PHYSICIAN: Ethylene oxide is a liquefied gas. Skin exposure from contact with fabric, rubber, or plastic containing residual ethylene oxide will commonly result in skin irritation with extensive blister formation. At high concentrations severe conjunctivitis can occur. Irritation of the respiratory tract may occur, but without acute lung edema. Symptoms of systemic intoxication are headache, nausea, vomiting, lack of coordination, and cardiac irregularities. Treatment is symptomatic.

PERSONAL PROTECTIVE EQUIPMENT (PPE): A material that is chemical resistant to this product is butyl rubber. All handlers must wear at a minimum:

- Long-sleeved shirt and long pants.
- Shoes plus socks.

- Chemical-resistant gloves and
- when the ambient ETO concentration is 1 to 50 ppm, full-facepiece respirator with ETO approved canister, front or back mounted,
- when the ambient ETO concentration is 50 to 2,000 ppm, (1) positive-pressure supplied-air respirator equipped with full-facepiece, hood, or helmet; or (2) continuous-flow supplied-air respirator (positive-pressure) equipped with hood, helmet, or suit,
- when the ambient ETO concentration is >2,000 ppm or unknown (e.g., emergency situations), (1) positive-pressure self-contained breathing apparatus equipped with full-facepiece; or (2) positive-pressure full facepiece supplied-air respirator equipped with an auxiliary positive-pressure self-contained breathing apparatus.

When handlers could have eye or skin contact with ETO or ETO solutions, such as during maintenance and repair, vessel cleaning, or cleaning up spills, they must wear:

- Chemical-resistant attire, such as an apron, protective suit, or footwear that protects the area of the body that might contact ETO or ETO solutions, and
- Face-sealing goggles, a full-face shield, or a full-face respirator.

USER SAFETY REQUIREMENTS

When wearing respirators:

1. Follow the respirator manufacturer's user's instructions for changing canisters.
2. Respirators must be fit-tested and fit-checked using a program that conforms to OSHA's requirements (see 29CFR Part 1910.134).
3. Respirator users must be trained using a program that confirms to OSHA's requirements (see 29CFR Part 1910.134).
4. Respirator users must be examined by a qualified medical practitioner to ensure physical ability to safely wear the style of respirator to be worn. A qualified medical practitioner is a physician or other licensed health care professional (PLHCP) who will evaluate the ability of a worker to wear a respirator. The initial evaluation consists of a questionnaire that asks about medical conditions (such as a heart condition) that would be problematic for respirator use. If concerns are identified, then additional evaluations, such as a physical exam, might be necessary. The initial evaluation must be done before respirator use begins. It does not need to be repeated unless the health status or respirator use conditions change (see 29CFR Part 1910.134).

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Follow manufacturer's instructions for cleaning/maintaining PPE. If no such instructions for washables exist, use detergent and hot water. Keep and wash PPE separately from other laundry

ENVIRONMENTAL HAZARDS

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollution Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product into sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

Emergency Contact: 1-800-255-3924

PHYSICAL AND CHEMICAL HAZARDS

Ethylene oxide gas is extremely flammable. Do not use near flame, electrical sparks, hot surfaces, or allow sources of ignition near the sterilization area. Ground all equipment to prevent static sparks. Contents under pressure. Do not puncture or incinerate container. Exposure to temperatures above 130 °F (54°C) may cause bursting.

ODOR: Ether-like in high concentrations. Exposure to toxic levels may occur without warning or detection by the user.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Employers in facilities that use ETO must comply with all requirements for ETO use specified in 29 CFR 1910.1047.

This product may be used only in Andersen AN74, AN75, EOGas 3, and EOGas 4 sterilizers for commercial, hospital, medical, and veterinary sterilization.

In hospitals, healthcare facilities, and contract sterilization facilities treating medical equipment and supplies, sterilization/fumigation with ETO must be performed only in vacuum or gas tight chambers designed for use with ETO. After February 28, 2010, a single chamber process is

required for ETO treatment (sterilization and aeration are to occur in the same chamber) in hospitals and healthcare facilities.

During the sterilization cycle, the sterilization bag may contain hazardous concentrations of ETO. Do not open the sterilization bag until the sterilization and ventilation portions of the sterilization cycle are complete.

Safety and awareness training is required for all employees including office staff. Information and training must be provided to all employees in the facility at the time of initial assignment and annually thereafter. The safety training must include, at a minimum, the following information: 1. the most recent monitored ambient levels of ETO in the facility; 2. the potential health effects from the levels of ETO in the facility; 3. the emergency response plan and how to respond in an emergency; 4. the availability of the Material Safety Data Sheet and other materials related to the health hazards of exposure to ETO. To reduce ambient levels of ETO, lengthy facility aeration is encouraged. Achieving an ambient level of 0.25 ppm (measured as a daily time-weighted average) greatly reduces potential long-term risks to employees not directly involved in the ETO applications. Air monitoring should include the entire facility including office space, break areas, and loading/unloading areas. Employers in facilities that use ETO must comply with all the requirements for ETO use specified in 29 CFR 1910.1047.

Sterilization Process:

- Remove one sterilization bag from the dispenser box. Place wrapped devices to be sterilized inside the sterilization bag along with required accessories (Humidichip, chemical indicators, biological indicator).
- Remove one cartridge from the dispenser box. Remove the adhesive tape, then remove the trigger guard. Place the cartridge on top of the wrapped items near the open end of the bag where it is accessible for cartridge activation. **DO NOT ACTIVATE THE CARTRIDGE.**
- On sterilizers with sterilization bag purge, insert the purge probe into the sterilization bag. Gather the open end of the sterilization bag around the purge probe and wrap the Velcro strap around both the sterilization bag opening and the purge probe, tightening it to completely seal the sterilization bag. On sterilizers without a purge probe (EOGas 3), the bag must be sealed with a heat sealer after removing excess air from the sterilization bag.
- Place the sealed sterilization bag inside the sterilizer. On sterilizers with sterilization bag purge, connect the quick release connector on

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the purge probe to the purge probe tubing. In the AN75 and EOGas 4 sterilizers, include a biological indicator in the BI receptacle on the purge probe.

- In Anprolene and EOGas 4 sterilizers with a purge probe, evacuate excess air from the sterilization bag using the initial bag purge, then follow the cycle selection and cycle verification prompts.
- **ACTIVATE THE CARTRIDGE:** grasp the cartridge through the wall of the sealed sterilization bag and fully depress the trigger button on the side of the cartridge to release ETO into the sterilization bag.
- Close the door to the sterilizer and start the sterilization cycle.
- Always lock the door of the sterilizer to prevent unauthorized removal of devices during the sterilization and ventilation portions of the cycle.
- After the sterilization and ventilation portions of the cycle in sterilizers without a purge probe (EOGas 3), open the sterilization bag while it is still in the cabinet, then shut the door and continue aeration. This will improve aeration efficiency and minimize the chance of exposure to residual EO within the sterilization bag.
- After the sterile devices have been removed from the sterilization bag, discard the sterilization bag and empty cartridge according to disposal instructions. **DO NOT REUSE EMPTY BOX, EMPTY CARTRIDGES, OR LINER BAGS.**

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

Pesticide Storage: Store in a cool, well ventilated area away from heat and direct sunlight. Store in accordance with 29 CFR 1910.1047.

Pesticide Disposal: Pesticide wastes are toxic. Improper disposal of excess pesticides (unwanted or expired ampoules) is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact Andersen Sterilizers, Inc. (336-376-3000), or your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

Container Disposal: Do not puncture or incinerate unused ampoules. After aeration, dispose of empty ampoules and used sterilization bags in sanitary landfill.

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