



An ISO 13485 certified company

Andersen Products, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 USA
Phone: 336-376-3000

Re: Useful life of Andersen sterilizers

The useful life of all Andersen sterilizers is ten (10) years from the date of manufacture. The Lot/Serial Number on the back of the sterilizer indicates the year the sterilizer was manufactured. For example, the first two numbers of the Lot/Serial Number 110659 indicate that the sterilizer was manufactured in the year of 2011. Useful life refers to the time period during which Andersen Sterilizers will maintain a spare parts inventory and provide service to repair the sterilizer.